



QUBÉ

THE COOL BEAUTY TOOL

BEAUTY AND PAIN RELIEF THROUGH COLD



ATTENTION

Store in a protected place.

Do not use if broken or damaged, dispose in accordance with local regulations. In case of skin contact with the gel contained, rinse with water without rubbing.

WHAT IS QUBÉ?

It is a Class I medical device consisting of a coolable and reusable cube.

WHAT IS QUBÉ USED FOR?

Support for aesthetic therapy and cryotherapy

FOR WHOM IS IT SUITABLE?

MEDICAL BEAUTY TREATMENTS

- filler and laser treatments
- lip treatments (tattoos, fillers, bioregeneration)
- Biostimulation and biorevitalization
- Aesthetic and regenerative medical treatments

- Adjuvant in the reduction of Eyelid edema
- Adjuvant in the reduction of fatty bags on the eyelids

BEAUTY TREATMENTS

- depilation
- Facial cleansing treatments
- Aesthetic tattoos

- Improvement of rosacea
- Adjuvant in the reduction of Eyelid edema
- Adjuvant in reducing the signs of fatty bags on the eyelids

Can the cube be used for a long period of time?

The cube can be used on the instructions of a physician and as a cooling device for non-medical purposes.

When is it not recommended?

It is not indicated for the treatment of traumatized body parts when and/or complicated by ischemic lesions or vasculitis, or in areas where sensory perception is impaired due to peripheral nervous system damage. The cold carrier should not be placed in

direct contact with open wounds, ulcers and surgical wounds as healing may be impaired.

TYPE OF APPLICATION

Put the product in the enclosed bag and put it in the freezer for at least 6 hours (the product absorbs the maximum amount of cold after 24 hours), remove from the freezer and

take the product out of the bag, apply to the affected area until the temperature of the cube corresponds to the skin or for the specified period of time.

HOW TO CLEAN

It can be cleaned with normal non-aggressive detergents and water, for disinfection use a solution of hydrogen peroxide, do not put in the oven or microwave.



SKU Q2202
LOT 2113230568
MFD 01/2022

DISTRIBUTED by
koome Care & Cosmetics SL
Carrer Zurbaran 1
Roses (Spain)

MANUFACTURED by
IRIANS Monaco
7 Av. Saint Roman
98000 Monaco
Principality of Monaco

